

## TECHNICAL DATASHEET TP05



### PNEUMATIC COMPRESSION DEVICE

<b>Medical Device</b>	class IIA	Annex IX rules 9 93/42 CEE
<b>Medical Device</b>	class II	Schedule I rules 9(1) SOR/98-282
<b>Medical Device</b>	class II	21 CFR part 870 § 870.5800



## TECHNICAL SPECIFICATIONS

<b>Weight</b>	6,5 kg	<b>Sizes</b>	40   31   11 cm
<b>Power</b>	44 W	<b>Voltage / Frequency</b>	220 V / 50 Hz
<b>Sound level</b>	49 dB<l<52 dB	<b>Pressure Range</b>	20 to 200 mB (15 - 150 mmHg)
<b>Working Time</b> (on hold)	0 to 30s	<b>Rest Time</b> (on hold)	0 to 60s
<b>Treatment Duration</b>	0 to 60 min	<b>Operating Temperature</b>	+0°C to +40°C
<b>Operating Atmospheric Pressure</b>	70kPa to 106kPa	<b>Operating Relative Humidity</b>	30% to 80%
<b>Storage Temperature</b>	-5°C to +45°C	<b>Storage Relative Humidity</b>	30% to 80%
<b>Storage Atmospheric Pressure</b>	70kPa to 106kPa		

IEC 60601-1 : 94,2 00547,02

Class II Device Type BF IP20

### Pressure Regulation

Automated Double Gradient Pressure  
(A.D.G.P. – International Patent)

### Programs range

Unique protocole with manual adjustment of Pressure, Working and Rest times and Session duration.

### Garments inflation

Two rows of 5 connectors allow simultaneous inflation of 1 or 2 accessories.



Accessories Components	Connection
Outside: Nylon / Inside: Polyurethane Closing method : Thermo-welding	Tubes on Garments: Fixed to avoid air leakage Tubes on Pump: Individual through numbered connectors
Garments list	Innovation
5 chambers Sleeve 5 chambers Boot 5 chambers Waist Boot 5 chambers Abdominal Strap	Easy-to-use Pressotherapy Device with fixed protocole and direct adjustments. Robust and reliable, it is Eureduc iconic product.

 Eureduc

19, rue de Clairefontaine - 78120 RAMBOUILLET - FRANCE

Tel. : +33 (0)1 34 61 01 61

info@eureduc.eu - www.eureduc.eu

## ACCESSORIES

Ref.	Description	Sizes	
M5PGM	5 Chambers Sleeve	Ø 58 cm x 77 cm	
B5PIM	5 Chambers Boot	Ø 80 cm x 85 cm	
H5P	5 Chambers Waist Boot	Ø 80 cm x 106 cm	
SPK	5 Chambers Abdominal Strap	L: 56 cm x l: 37 cm	
GHB	Hygienic sheathes (100 units)	40 x 120 cm	
PU04	Garments Support	H: 60 cm	
PU09	Mobil Stand with Garments Support	l=41 cm x P=41cm x H=79cm	

### Indications

- Œdemas: Lymphœdema, Phlebœdema,
- Prevention of Deep Venous Thrombosis (DVT)
- Post-thrombotic syndrome
- Chronic Venous Insufficiency (CVI) symptoms / Heavy legs
- Recovery after sports or surgery

### Contra-indications

- Acute DVT
- Acute infection: Erysipela or Lymphangitis
- Systemic decompensation
- Pressure urticaria