


TECHNICAL DATASHEET TP05 i



PNEUMATIC COMPRESSION DEVICE

Medical Device	class IIA	Annex IX rules 9 93/42 CEE
Medical Device	class II	Schedule I rules 9(1) SOR/98-282
Medical Device	class II	21 CFR part 870 § 870.5800



TECHNICAL SPECIFICATIONS

Weight	9 kg	Sizes	37 33 18 cm
Power	40 W	Voltage / Frequency	100 to 240 V / 47-63 Hz
Sound level	49 dB<l<52 dB	Pressure Range	20 to 150 mBar (15-112 mmHg)
Medical Power Supply	GSM90A24-PIM Meanwell	Treatment Duration	0 to 60 min
Working Time (on hold)	0 to 60s	Rest Time (on hold)	0 to 60s
Operating Temperature	+0°C to +40°C	Operating Relative Humidity	30% to 80%
Operating Atmospheric Pressure	70kPa to 106kPa	Storage Temperature	-5°C to +45°C
Storage Relative Humidity	30% to 80%	Storage Atmospheric Pressure	70kPa to 106kPa

IEC 60601-1 : 130879-663528

Class II Device Type BF IP20

Pressure Regulation

Automated Double Gradient Pressure
(A.D.G.P. – International Patent)

Programs range

9 presets protocoles
+ 3 customized programs

Garments inflation

Two rows of 5 connectors allow simultaneous inflation of 1 or 2 accessories.



Accessories Components	Connection
Outside: Nylon / Inside: Polyurethane Closing method : Thermo-welding	Tubes on Garments: Fixed to avoid air leakage Tubes on Pump: Individual through numbered connectors
Garments list:	Innovation
5 chambers Sleeve 5 chambers Boot 5 chambers Waist Boot 5 chambers Abdominal Strap	Easy-to-use Pressotherapy Device with preset antegrade and retrograde protocoles controlled through touch screen

Eureduc

19, rue de Clairefontaine - 78120 RAMBOUILLET - FRANCE

Tel. : +33 (0)1 34 61 01 61 - Fax : +33 (0)1 34 61 40 61

info@eureduc.eu - www.eureduc.eu

ACCESSORIES

Ref.	Description	Sizes	
M5PGM	5 Chambers Sleeve	Ø 58 cm x 77 cm	
B5PIM	5 Chambers Boot	Ø 80 cm x 85 cm	
H5P	5 Chambers Waist Boot	Ø 80 cm x 106 cm	
SPK	5 Chambers Abdominal Strap	L: 56 cm x l: 37 cm	
GHB	Hygienic sheathes (100 units)	40 x 120 cm	
PU04	Garments Support	H: 60 cm	
PU09	Mobil Stand with Garments Support	l=41 cm x P=41cm x H=79cm	

Indications

- Œdemas: Lymphœdema, Phlebœdema,
- Prevention of Deep Venous Thrombosis (DVT)
- Post-thrombotic syndrome
- Chronic Venous Insufficiency (CVI) symptoms / Heavy legs
- Recovery after sports or surgery

Contra-indications

- Acute DVT
- Acute infection: Erysipela or Lymphangitis
- Systemic decompensation
- Pressure urticaria